

Wirral & West Cheshire Wheelchair Service

Correct Procedure for Connecting and Disconnecting Powered Wheelchair Battery Charger

Connecting the battery charger: -

- Switch off wheelchair
- Plug battery charger into the socket on the wheelchair
- Plug the battery charger into the mains socket and switch on
- Leave on charge overnight

Disconnecting the battery charger: -

- Check that the fully charged indicator light is lit on the battery charger. If it is not, contact the repairers and report the fault
- Switch the battery charger off at the mains
- Disconnect the battery charger from the wheelchair

Notes

1. Connecting or disconnecting the battery charger from the wheelchair whilst the mains are still on will, over a period of time, result in damaged connections within the wheelchair battery charging socket. This will result in the wheelchair having a reduced range and reduced power.
2. This information leaflet does not replace the information provided by the manufacturer with the wheelchair.