


WIRRAL & WEST CHESHIRE WHEELCHAIR SERVICE**IMPORTANT NOTICE****Invacare Alulite Wheelchairs**

This wheelchair has not been designed or tested to be used as a seat in a vehicle.

This means you cannot sit in your wheelchair in ambulances or adapted vehicles.

However traveling in public transport such as buses or trains (where the wheelchair is not tied down) is permitted.

If your wheelchair needs to be tied down in a vehicle please contact the Wheelchair Service and we can reassess you for an alternative wheelchair.